











How to set up your replacement MiniMed™ 780G

These videos will help you to enter your settings, and with additional features

Before you start, we recommend that you download your current settings from your pump, by logging onto your [CareLink™ Personal](#) account. For more information regarding accessing your CareLink™ Settings Report, please click [here](#). Alternatively you can use our [Settings Sheet](#) to record your settings. Please scan the QR code on the right to download this.

Use the Startup Wizard to set the language, time format, current time and date, and to rewind the pump, and then click on each video to view the step by step guides for entering your settings into your new pump.



Deleting Paired Devices 	Connecting BG Meter 	Connecting Transmitter 	
Setting Basal Pattern(s) 	Setting up Bolus Wizard 	Setting Max Basal and Max Bolus 	Sensor Settings 
Entering a BG 	Suspending and Resuming Insulin 	Using Temp Target 	Using Temp Basal 